

# VETSCAN Rotor Guarantee Report

Vet Practice Address: \_\_\_\_\_

Eircode: \_\_\_\_\_ Customer Contact: \_\_\_\_\_

Practice E-mail address: \_\_\_\_\_ Phone: \_\_\_\_\_

VETSCAN VS2 Analyzer: SN Serial Number: \_\_\_\_\_

SW version: \_\_\_\_\_

Profile Types (please tick the failed one) + Lot Number: \_\_\_\_\_

- ☐ Phenobarbital Profile
- ☐ Electrolyte Plus
- ☐ Comprehensive Diagnostic Profile
- ☐ Prep Profile II
- ☐ Mammalian Liver Profile
- ☐ T4/ Cholesterol Profile
- ☐ Equine Profile Plus
- ☐ Kidney Profile Plus
- ☐ Avian/ Reptilian Profile Plus
- ☐ Large Animal Profile

Sample Type: ☐ Whole Blood ☐ Serum/Plasma

Sample Collection Tube: ☐ Lithium Heparin ☐ Other \_\_\_\_\_

iQC Code: \_\_\_\_\_

Date of Submission: \_\_\_\_\_

Additional  
Comments: \_\_\_\_\_

Thanks for your time in completing this form!

Please send this filled in form to [vetscanireland@zoetis.com](mailto:vetscanireland@zoetis.com)

Please contact the Technical Support on 0044 845 300 8034 for any emergency technical support.

## Please stick your results here

(including the iQC report)  
(You may remove any patient data however all other  
information should be clearly visible please)

The VETSCAN VS2 has an intelligent Quality Control  
(iQC) that automatically verifies the integrity of samples,  
chemistries and analyzer performance  
during each run. It automatically suppresses either  
a single chemistry or the entire panel if it detects any  
endogenous substance interference by printing flags  
instead of results  
(~~~, HEM, LIP or ICT, iQC Code)

Moreover, the run can be cancelled due to a rotor failure.

This ensures that accurate results are reported,  
and it prevents misinterpreting of  
results that could be affected  
by sample quality issue

In line with GDPR guidelines we will not process any  
patient data, please follow your own practice GDPR  
policies when sharing data

[All trademarks are property of ZOETIS Services LLC or a related company or licensor unless otherwise noted . © 2019 ZOETIS LLC . All rights reserved](#)

MM-09646

**zoetis**